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NOTICE: This document is version controlled and was produced as a part of the GEX Information Program which requires that all Series 100 documents be reviewed periodically to maintain currency and continuity of information. Appropriate Technical Memorandum are used to provide information detail in support of the Product Data Sheets as well as GEX Recommended Procedures and to provide technical information in support of GEX Marketing documents.

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This procedure describes the methods of evaluating the uncertainty of instrumentation and personnel during dosimeter calibration and routine use.

2.0 材料 MATERIALS

2.1 WINDose 剂量测量系统
WINDose Dosimetry System

2.2 GEX Doc#100-252, 仪器和人员鉴定表
GEX Doc#100-252, Instruments and Personnel Characterization Form

3.0 频率 FREQUENCY

3.1 一批剂量计剂量响应鉴定期间
During the dose response characterization of a batch of dosimeters

3.2 培训新员工
Training of new hires

3.3 从业人员年度鉴定
Annual re-certification of employees

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Remove the cuvette cup from the sample compartment and close the Genesys 20 Spectrophotometer sample compartment cover. Turn the instrument on from an off position and allow the instrument to go through its complete automated start-up sequence that will take about 3 minutes to complete. After the instrument has completed its start-up sequence, allow a full 30 minute warm-up period before attempting to use the instrument. Following completion of the warm-up period it is safe to replace the cuvette cup and zero the instrument at the appropriate absorbance measurement wavelength for use.

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If using the WINdose dosimeter holder, insert the holder with the conical facing toward the light beam and zero the instrument. Practice inserting the holder several times and achieving “zero” readings until the comfortable with the placement technique.

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Begin documenting the individual operator’s results. Insert the empty holder 32 times and record the readings (including positive and negative absorbance values) on the Instruments and Personnel Characterization Form.

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Successful performance of the test is met when no measurement greater than $\pm 0.003A$ or greater is observed during the 32 cycle test.

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MEASUREMENT INSTRUMENTATION VARIANCE TESTING

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Select a dosimeter that has been irradiated to approximately the middle of the range 25 and 50 kGy (or the middle portion of the calibrated range).

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Measure the irradiated absorbance, A_i , of the test film by placing the dosimeter into the film holder positioning the holder into the sample cell holder of the Genesys 20, closing the lid and recording the displayed absorbance value. Open the sample compartment lid, remove the dosimeter holder from the instrument and repeat this cycle a minimum of 32 times in each spectrophotometer. Enter these values into the appropriate cells in the Instruments and Personnel Characterization Form or equivalent.

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Remove and replace the dosimeter from the holder each time. This tests typical handling uncertainty during routine dosimeter measurements. Use care when handling the dosimeter so as to avoid dents, smudges, or the accumulation of particulate on the film during this test. A damaged or dirty dosimeter film may skew the test result.

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The preformatted formulae in the Instruments and Personnel Characterization Form spreadsheet will calculate the mean, minimum, maximum, standard deviation, and coefficient of variance (CV) for the readings from the test films. Performing three consecutive tests for a single operator following the instruction in 5.1, 5.2 and 5.3 provides data necessary to determine an estimate of the measurement instrumentation (spectrophotometer and dosimeter holder) repeatability (single dosimeter, single instrument system and single operator). This test can be repeated over several days with different operators to establish a reproducibility **variance** for the instrument system. The data collected from this testing can be used to establish a baseline for the instrument system. This information is also useful in understanding the contribution of the dosimeter system to overall dosimetry system uncertainty.

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The instrument system should be expected to result in a repeatability CV of 0.3% or less at an absorbance of 0.250 or higher if using a DoseStix holder and 0.5% when using the WINdose holder.

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The dosimeter may be discarded upon completion of testing.

6.0 修订历史 REVISION HISTORY

日期 Date	版本 Revision	修改描述 Change Description
07/27/07	C	增加了明确的仪器预热说明 Added explicit instrument warm –up instruction 修改了关于系统差异的语言描述 Changed language pertaining to system variance 修改了对于系统差异的测试通过/失败标准 Changed test pass/fail criteria for system variance.